

movinglife®

are created specifically to enhance, support and work with the ATTO mobility scooter in all user modes, while driving, while wheeling along as trolley and when folded for storage or transportation.

**ATTO ACCESSORIES**

**ATTO**  
MORE FOR

**ATTO**  
CARRYALL & CUSHION


Don't Think Twice,  
**JUST BRING IT ALONG!**

movinglife®



## HOW TO USE

# ATTO CARRYALL & CUSHION



### 01


Slip the cushion over the seat (use handles to pull if needed)



### 02

Simply pull hook on elastic strap to fold under and connect to ATTO chassis

**\* NOTE!** Take off before folding ATTO



### 03

One-hand "pull to open" secure funnel access, so nothing falls out.

Quick "no zip" flap, to ensure everything stays inside when closing.